Bruksanvisning Operating Instructions

Master-care[®]
The Swedish Back Care System[™]

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
This instruction manual should always be kept available by the
Mastercare Back-A-Traction Mini.

Please read the manual carefully and follow the instructions.

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Before you can start using Mastercare Mini you should remove the transit bolt. This bolt should then be located in the vertical slot situated directly above it. It should be aligned with the matching hole situated on the sliding frame beneath the backboard. There are four vertical slots on the sliding frame and this hole is situated at the top of the bottom slot.

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
Pictures and descriptions in this publication were correct at time of printing.

We reserve the right to make changes due to a continuous development and improvement of our products.

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ENGLISH

Back-A-Traction Mini is a gravity traction table that will help you relieve and prevent back pain and the wearing out of joints. It is a good tool to help you regain your back health, and to maintain it for the rest of your life.

Back-A-Traction Mini's unique feature the Sliding Backrest makes it possible to obtain optimum beneficial traction at a low inversion angle. While your spine is extended, the back is perfectly supported allowing the muscles to relax completely. You will experience an unloading of the pressure of gravitational forces from your discs and joints.

To get optimum results, you should alternate stretching exercises and total relaxation. The passive traction and the stretching exercises on the table involve stretching of the muscles, muscular attachments, tendons and crucial ligaments of your entire back and legs. In the beginning this may cause a slight stiffness, which is a sign that the mentioned parts of your body are not used to this kind of stretching.

Even if you become sore in the beginning, you may continue your exercises with less intensity, as the soreness should disappear as the muscles and tendons increase their flexibility. Daily use of Back-A-Traction Mini will lead to a healthier life.

CAUTION ADVISED IN CERTAIN MEDICAL PRE-EXISTING CONDITIONS

ALWAYS consult your doctor before starting a regular use of the table if you are uncertain about the state of your health, or suffer from any of the following conditions:

Glaucoma, detached retina, conjunctivitis, uncontrolled hypertension (high blood pressure), impaired pulmonary functions, arrhythmia, carotid artery stenosis, severe vascular disease, aneurysm, tendency towards fainting or dizzy spells, cerebral sclerosis, recent stroke or transient ischemic attack, history of brain lesions, other heart or circulatory disorders for which you are being treated, severe back problems or spinal injuries, bone weakness (osteoporosis), recent unhealed fractures, medullary pins, artificial hip joints or other orthopedic supports surgically implanted, middle ear infection, recently had surgery, using special medication, ventral or hiatal hernia, tendency towards bleeding, extreme obesity, invertophobia (fear of inversion itself), and pregnancy.

IMPORTANT SAFEGUARDS

Not for use by minors themselves under the age of 18.

Have others remain clear while the table is in use, to avoid bodily injury from moving parts.

As with any therapy or exercise equipment, you should consult your doctor before regular use of Back-A-Traction Mini.

The table is designed for use by one person at a time, with both feet fastened in the foot-locking device. In the beginning, use Back-A-Traction Mini in the presence of a person who is familiar with its functions.

Before using Back-A-Traction Mini on your own, you must be certain that you fully understand how the table works.

Be sensitive to the signals from your body. Consult your doctor or physical therapist if you feel uncomfortable, have pain or notice any unusual changes that you think might prevent you from using Back-A-Traction Mini.

Do not use Back-A-Traction Mini immediately after a meal, or if you have a temperature, a cold or severe head- or ears ache.

Back-A-Traction Mini is designed for persons with a height between 4,6 and seven feet (135 – 215 cm), and a maximum weight of 220 lbs (100 kilos).

Make sure that:

- it stands on a level, horizontal surface
- all nuts and bolts are securely fastened
- the foot lock is working perfectly

Learn how the thumb catch works and make sure that it is engaged and working properly.

Before fastening your feet in the foot holder, always make sure that the catch for the height adjustment is in a locked position.

If you do not balance well when operating Back-A-Traction Mini, re-adjust the foot rest to obtain the best possible balance. (In order to make sure that the backrest will not rock back and forth when you do the stretching exercises, it is usually advantageous to have it adjusted so that the part where your upper body is resting weights a little more). Wiggle the backrest slightly when you loosen the safety catch. If the thumb catch works stiffly, the height is not properly adjusted. To correct this, the foot rest has to be extended.

If your back is in a very bad condition when you start using Back-A-Traction Mini, it will for the first couple of weeks, probably be enough for you to just hang in it, with your arms over your head. Then, as your condition improves, you can begin to do the stretching exercise program.

A BASIC STRETCHING EXERCISE PROGRAM

Follow the instructions you have received. Be cautious and sensitive to the signals from your body. For best result you should relax and take some deep breaths between the different stretching exercises. It is also important that you relax as much as possible during each stretching exercise.

The following instruction will tell you to repeat certain exercises. It is for you to decide if you want to repeat them, and how many times. An easy start is recommended. Increase muscular stretching and relaxation with each repetition. Most of the following stretching exercises should be done, if possible, with your arms over your head.

1. Start with unfolding the stand by pressing the big U-shaped plastic button and pushing the rear end of the stand (picture 2 and 3). After that it is time to adjust the balance of the back rest according to your height.
2. Press the thumb lock and rotate the back rest to the horizontal position where it locks in place (pictures 4 and 5).
3. To adjust the balance, pull the height and balance adjustment knob and slide the bar according to the guideline markings provided. Test the back rest and fine tune the balance to your height and weight as necessary (pictures 6 and 7). Then adjust the gap in the foot holder to fit your feet (picture 8).
4. Press the thumb lock and rotate the back rest towards its vertical position. It will lock at a slight incline. (picture 9).

You're now ready to use Back-A-Traction Mini. As a first time user you should always have someone help you until you feel comfortable using Back-A-Traction Mini.

5. Remove your shoes. Sit on the lower end of the back rest while placing your feet between the foot holder pads. Check that your feet will stay and are held, comfortably and securely between the pads (pictures 10 – 12).
6. Now, lean back onto the back rest and take hold of the hand grips and with your right hand thumb on the button and press. Release the thumb lock and push back to the horizontal position where the back rest will lock in place (pictures 13 and 14).

Before exercising, you should stretch by pushing against the hand grips (picture 15). Throughout the session, remember to relax and breath deeply and allow Back-A-Traction® Mini™ to work, stretching the tension from your muscles, tendons and ligaments, the common causes of back pain. (You may adjust the “moving space” of the back rest by changing your position on the board).

7. After stretching, press and release the thumb lock and rotate further backwards to the catch at 15 degrees. Repeat the stretching procedure by pushing, as before, against the hand grips. Now, raise your hands slowly over your head and prepare to exercise. First begin with a few repetitions of bending your feet under the restraint pads (pictures 16 and 17).
8. Next, alternate bending your knees in a slow and easy walking motion. Each exercise should be repeated a few times for best effect (pictures 18 and 19).

Remember to relax and breathe deeply between each exercise.

9. The exercise involves arching and lifting your buttocks off of the back rest (picture 20).
10. Next, push your neck flat to the back rest while lifting your shoulders upward (picture 21).
11. After a few repetitions, place one hand on your stomach and then reach over your head with the other, as if trying to pick apples just out of reach (picture 22).

Now, repeat the process, this time holding your arm out to your side. You should alternate your palms up and down from time to time during this exercise. Finally, repeat the exercise holding your hand near your waist (picture 23).

From this finishing position, move right into the next exercise, moving your arm in circles as if swimming the backstroke. Repeat these hand and arm exercises with your other arm. End with your back flat to the backrest and both hands together on your stomach.

12. The next exercise involves a three part move. Raise both hands over your head and then twist your torso and hips to the right, raising the right knee at the same time. Return to the flat back, flat knee position and twist to the left side. Return to the flat position and repeat the move a few times. End the exercise by bringing your hands back to rest on your stomach (picture 24).
13. The final exercise involves quarter sit-ups. Raise your head, shoulders, and upper back only, leaving the rest of your back flat to the backrest. Proceed only as you feel comfortable to do so (pictures 25 and 26).
14. Clasp your hands together on your stomach. Turn your head as far as possible towards each side. Relax. Take some deep breaths. Repeat (picture 27).

End by returning flat to the backrest and raising your arms over your head for a last reaching stretch.

15. Begin the dismount by taking hold of the hand grips, releasing the thumb lock, and returning to the horizontal position. Relax for at least 60 seconds and move your knees in a walking mode for a few repetitions as before.
16. Finally, release the thumb lock and return the back rest to its last and locking point. Release your feet from the foot pads (pictures 28 and 29).

Step out and away from your Back-A-Traction Mini (picture 30).

It is very important to relax after each session. Avoid lifting and strenuous movements, and take it easy for a few minutes.

Now you have a basic stretching exercise program.

Start with caution. In the beginning "hold" means a few seconds. After you have been training for some weeks it may mean 20 to 30 seconds. Notice the reactions of your body. Intensify the exercises gradually.

Feel free to talk to your back specialist about an exercise program just for you.

Total relaxation is necessary to get the best results.

The most important thing you do on the table is to relax, to let your muscles and tendons stretch. The muscles let "go" of their tension little by little as you "hang" and relax, and by this you regain your flexibility. There is also a beneficial "negative" pressure on your discs, a "suction" action, achieved from just "hanging".

We recommend that you use Back-A-Traction Mini one to three times per day, two to fifteen minutes per session.

It is important that you relax and take it easy for a few minutes after each session. You may find it very beneficial just to hang and relax during an entire session. We also want you to be prepared that you may experience "popping" in your back as your spine is being aligned.

Please visit www.mastercare.se for more information.